



Radioactive Contamination

Monitoring Methods



Two Methods

- Hand-held monitor
 - Usually a Mini-Instruments mini-monitor
 - Geiger-Muller and solid scintillation detectors
- Wipe test
 - Liquid scintillation counting

Standard Mini-Monitor

Battery compartment

Off/Bat/On/Mute switch

Alarm level set

Battery indicator

Rate-meter scale


Probe

End window



Type E/S Probes

- Primarily used to detect β radiation
- Geiger-Muller (GM) detectors






Type E/S Probes

E	<p>6cm² end-window. Background = ~0.25cps. Used for hard β, e.g. ³²P, ³⁶Cl. Will detect soft β, e.g. ¹⁴C, ³⁵S, ³³P, ⁴⁵Ca, but EP15 is preferred. Will not detect very weak β, e.g. ³H. Sensitive to γ, X and α.</p>
EP15	<p>Twice as sensitive as type E. 15cm² end-window. Background = ~0.5cps Used for hard β, e.g. ³²P, ³⁶Cl. Especially good for soft β, e.g. ¹⁴C, ³⁵S, ³³P, ⁴⁵Ca, except ³H. Also sensitive to γ and X. OK for $\alpha > 3$Mev.</p>
S/SL	<p>Open grille on one side. Background = ~1cps. Used for hard β, e.g. ³²P, ³⁶Cl. Not suitable for soft β, e.g. ³H, ¹⁴C, ³⁵S, ³³P, ⁴⁵Ca. Suitable for β/γ emitters such as ²²Na.</p>



^{22}Na Decay Scheme





β Energies

Isotope	βE_{\max}
3H	0.019
14C	0.156
35S	0.167
33P	0.249
45Ca	0.257
36Cl	0.710
32P	1.711

Not detectable by mini-monitor

Soft β


EP15

Hard β

E or S/SL

Type 42/44 Probes

- Used to detect γ and x-ray radiation
- All NaI solid scintillation detectors





Type 42/44 Probes

42A	0.05mm thick aluminium end window. 23mm diam, 1mm thick crystal. Background = 1.5 to 3 cps. Used for low intensity γ/x , 10 to 150 KeV.
42B	0.25mm thick beryllium end window. Much better detection efficiency at energies < 10 KeV. Used for low intensity γ/x , 5 to 150 KeV.
44A	0.05mm thick aluminium end window. 32mm diam, 2.5mm thick crystal. Background higher than Type 42 but much better efficiency above 40KeV. Used for wider range of γ/x , 10 to 500+ KeV.
44B	0.25mm thick beryllium end window. Much better detection efficiency at energies < 10 KeV. Used for wider range of γ/x , 5 to 500+ KeV.




Detector Efficiency

- % Efficiency = $\text{cps/dps} * 100$
- Determined by
 - Geometry
 - Absorption in air and end window
 - Intrinsic efficiency of the detector

Bremsstrahlung

Braking radiation





Using a Mini-Monitor

- Check
 - Next test date
 - Probe type
 - Battery strength
- Remove end-window cap
- Check
 - Speaker is on
 - Background reading
 - Response to check source



Monitoring Protocol

- Note date and type of monitor/probe used
- Note background cps
- Monitor slowly and methodically
- Hold the probe $\sim 1\text{cm}$ above the surface
 - Take care not to contaminate the end of the probe!
- Record
 - cps if $>$ background
 - zero if cps = background



Dealing with Contamination

- Establish the extent & demarcate area
- Attempt to decontaminate
- Continue until all loose contamination has been removed
- Re-monitor area
- If bound contamination remains and presents a significant external radiation hazard, shield it



Wipe Testing

- More laborious than mini-monitoring
- Essential for tritium
- Employs liquid scintillation counting
- Low efficiency
 - Usually only 10% of loose contamination is picked up
 - Subject to quenching
- Qualitative rather than quantitative



Possible Protocol


- Label counting vials (on cap not side)
- Use EtOH soaked filter disc, held in forceps
- Swab 10 x 10cm area (if possible)
- Work from the perimeter inwards
- Transfer disc to counting vial
- Take background wipe from an area never used for radioactive work, e.g. office bench
- Add 10ml of scintillation cocktail to each vial
- Shake well
- Count for 5 min on scintillation counter
 - Parameters will depend upon type of counter
- Any result 2x background is positive
 - Decontaminate and re-test until no loose contamination is detected



Alternative Methods

- Use a cotton-bud for very small areas
 - Break off into a microfuge tube
 - Add 150 μ l of scint fluid
- Use 2.5 diameter filter disc
 - Insert vertically into minivial
 - Add 0.5ml of scint fluid
 - Tilt carefully to completely wet filter
 - Insert into counter so disc faces detector

Liquid Scintillation Counting





Quenching

- Photon – β is absorbed
- Chemical – blocks solvent/fluor transfer
- Optical
 - Inside vial - colour quenching
 - Residue on outside of vial
- Other phenomena that interfere with LSC
 - Chemiluminescence
 - Static electricity